

CCK-8 Cell Counting Kit

Catalog # A311



Version 8.1

Vazyme biotech co., ltd.

1. Introduction

CCK-8 Cell Counting Kit is based on WST-8 [2-(2-methoxy-4-nitrophenyl)-3-(4-nitrophenyl)-5-(2,4-disulfophenyl)-2H-tetrazolium, monosodium salt] that is rapidly absorbing, highly sensitive and widely used in detection of cell proliferation and cytotoxicity. WST-8, similar to MTT, reduces to a water-soluble formazan dye in the presence of the electron carrying dehydrogenases found in mitochondria. The amount of formazan dye generated by dehydrogenases in living cells is directly proportional to the number of living cells. Measure the Optical Density at 450 nm, using a microplate reader. The measurements can be indirectly shows the number of living cells. CCK-8 Cell counting Kit is a ready-to-use solution, can be directly added to cell supernatants, and incubated for a certain period of time and then test.

2. Contents of Kit

Component	A311-01 (500 rxn, 10 µl / rxn)	A311-02 (1000 rxn, 10 µl / rxn)
CCK-8 Solution	5 ml	10 ml

3. Storage

Store at 4°C and protect from light.

4. Notes

1. CCK-8 Cell Counting Kit is pink solution, please protect from light.
2. The incubation time varies by the type and number of cells in a well. Please set several parallel well to make sure the appropriate amount of cell and the incubate time of CCK-8 solution. (Generally, the incubate time should be 1h to 4h)
3. Be careful not to introduce bubbles to the wells, since they interfere with the O.D. reading.
4. Since the CCK-8 assay is based on the dehydrogenase activity detection in viable cells, reducing agents (such as antioxidant) interfere with testing. Please try to remove reductant before using CCK-8 Solution.
5. In the experiments of pharmacological inhibition, if the metal is contained in drugs, such as Pb^{2+} , Fe^{2+} , Cu^{2+} , will interfere the color reaction of CCK-8 Solution, leading to lower the sensitivity of detection.
6. When measuring cell number, in order to ensure the stability and repeatability of the test results, it recommended drawing a standard curve at the same time.

5. Required Equipments and Materials

100 - 200 µl multi-channel pipettes, plate reader (with a 450 nm filter), 96-well plates, CO₂ incubator.

6. Protocol

This kit can be used for cell proliferation induced by cytokines and cytotoxicity test caused anti-cancer drugs, and drug-induced cell growth inhibition test.

1. Standard Curves

- 1) Collect and culture viable cells, calculate the number of cells in cell suspension with blood count plate, and then seed cells.
- 2) Inoculate cell suspension (100 µl / well) in a 96-well plate, dilute cells with medium in equal ratio gradient (such as 1:2). Five to 7 cell concentration gradients are recommended, with 3 to 6 replications in each group. Add 100 µl of cell suspension to each well.
- 3) Add 10 µl of the CCK-8 solution to each well of the plate. Incubate the plate in the incubator for an appropriate length of time. Measure the absorbance at 450 nm using a microplate reader. Take cell number as the abscissa, absorbance values as the ordinate to draw standard curve. Cell number of unknown samples can be measured according to the standard curve. The precondition of using the standard curve is that the test conditions is exactly the same.

2. Cell Number Determination

- 1) Inoculate cell suspension (100 µl / well) in a 96-well plate. Pre-incubate the plate in a humidified incubator (e.g., at 37°C, 5% CO₂) for 24 h. Set the blank group and the control group at the same time.
- 2) Add 10 µl of the CCK-8 solution to each well of the plate. (Be careful not to introduce bubbles to the wells, since they interfere with the O.D. reading).
- 3) Incubate the plate for 1 - 4 hours in the incubator.
- 4) Measure the absorbance at 450 nm using a microplate reader.



Vazyme Biotech Co., Ltd
www.vazyme.com

Order: global@vazyme.com

Support: global@vazyme.com

For research use only, not for use in diagnostic procedures.

3. Cell Proliferation and Cytotoxicity Assay

- 4) Dispense 100 μ l of cell suspension (5000 cells/well) in a 96-well plate. Pre-incubate the plate for 24 hours in a humidified incubator (e.g., at 37°C, 5% CO₂). Set the blank group and the control group at the same time.
- 5) Add 10 μ l of various concentrations of substances to be tested to the plate. Incubate the plate for an appropriate length of time (e.g., 6, 12, 24 or 48 hours) in the incubator.
- 6) Add 10 μ l of CCK-8 solution to each well of the plate. (Be careful not to introduce bubbles to the wells, since they interfere with the O.D. reading.) Incubate the plate in the incubator for an appropriate length of time.
- 7) Measure the absorbance at 450 nm using a microplate reader.

4. Calculating Formula

Cell viability % = $[(A - C) / (B - C)] \times 100\%$

Inhibition % = $[(B - A) / (B - C)] \times 100\%$

A: Experimental group OD (contain medium, cells, drugs and CCK-8 solution)

B: Control OD (contain medium, cells, CCK-8 solution)

C: Blank OD (contain medium, CCK-8 solution)

