

Product Description

AceTaq DNA Polymerase is a chemically modified Taq DNA Polymerase that is completely blocked at room temperature and is released only after heating at 95°C. Non-specific amplification and primer dimerization can be prevented during sample preparation and temperature rise. Compared with the antibody-based hot-start Taq, the polymerase activity of AceTaq DNA Polymerase is blocked more stringently and completely. It takes only 5 min to activate AceTaq DNA Polymerase. AceTaq DNA Polymerase is compatible with most existing PCR protocols. Combined with an optimized buffer system, AceTaq DNA Polymerase minimizes non-specific amplification and primer dimers, ensuring extremely high sensitivity and specificity, which make it ideal for amplifying low-copy genes from complex templates. The PCR products contain A at the 3' end and can be directly cloned into T-Vector. The products are compatible with ClonExpress and Top Cloning kits (Vazyme #C112/C113/C115/C601).

Components

Components	P401-d1	P401-d2	P401-d3
10 × AceTaq Buffer (Mg ²⁺ Plus)	1 ml	4 × 1 ml	
dNTP Mix (10 mM each)	200 µl	800 µl	3 × P401-d2
AceTaq DNA Polymerase (5 U/µl)	50 µl	200 µl	

Storage

Store at -30 ~ -15°C and transport at ≤0°C.

Applications

This product is suitable for DNA amplification from various type of templates such as genomic DNA, cDNA, plasmid DNA and λDNA.

Source

It is cloned from *Thermus aquaticus* and purified from *E. coli*.

Unit Definition

One unit (U) is defined as the amount of enzyme that incorporates 10 nmol of dNTP into acid-insoluble material in 30 min at 74°C with activated salmon sperm DNA as the template/primer.

Notes

For research use only. Not for use in diagnostic procedures.

Experiment Process

Reaction system

ddH ₂ O	to 50 µl
10 × AceTaq Buffer (Mg ²⁺ plus)	5 µl
dNTP Mix (10 mM each)	1 µl
Template DNA ^a	x µl
Primer 1 (10 µM)	2 µl
Primer 2 (10 µM)	2 µl
AceTaq DNA Polymerase (5 U/µl) ^b	0.5 µl

a. The optimal concentration for various templates is different. The recommended amount of DNA template for a 50 µl reaction is as follows:

Human genomic DNA	1 - 500 ng
<i>E. coli</i> genomic DNA	1 - 100 ng
λDNA	0.1 - 1 ng
Plasmid DNA	0.1 - 1 ng

b. The amount of AceTaq DNA Polymerase can be adjusted between 0.25 - 1 µl. Increasing the amount of enzyme may increase the yield, but it may also reduce the specificity.

PCR Program

95°C	5 min ^a (initial-denaturation)	
95°C	30 sec	} 30 - 35 cycles
55°C ^b	30 sec	
72°C	60 sec/kb	
72°C	7 min (final extension)	

a. The initial denaturation takes at least 5 min. If the amplification is not ideal, extend the initial-denaturation time up to 10 min.

b. Annealing temperature is based on the T_m value of the primers and is generally 3 ~ 5°C lower than the calculated T_m value.

Primer Design Guidance

1. It is recommend that the last base at the 3' end of primer should be G or C.
2. Consecutive mismatches should be avoided in the last 8 bases at the 3' end of the primer.
3. Avoid hairpin structures at the 3' end of the primer.
4. Differences in the T_m value of the forward primer and the reverse primer should be no more than 1°C and the T_m value should be adjusted to 55°C to 65°C (Primer Premier 5 is recommended to calculate the T_m value).
5. Extra additional primer sequences that are not matched with the template, should not be included when calculating the primer T_m value.
6. Control the GC content of the primer to be 40% - 60%.
7. The overall distribution of A, G, C and T in the primer should be as even as possible. Avoid using regions with high GC or AT contents.
8. Avoid the presence of complementary sequences of 5 or more bases either within the primer or between two primers and avoid the presence of complementary sequences of 3 or more bases at the 3' end of two primers.
9. Use the NCBI BLAST function to check the specificity of the primer to prevent non-specific amplification.

*All trademarks are the property of their respective trademark owners. Some trademarks are not registered in all administrative regions.