

1 × Phosphate Buffered Saline (PBS)


G101

Version 21.1



Product Description

On the basis of the combination of biology and pharmacy, Instant premixed granules series products have realized the large-scale automatic production process through technological innovation and production process transformation and upgrading, which greatly improved the stability of the product and ensures the repeatability of the experiment. Its application in conventional experiments can make the experimenter get rid of the shackles of balance, weighing paper and pH meter, which can greatly save time and carry out the experiment quickly and efficiently.



Components

10 mM Phosphate Buffer, 137 mM NaCl and 2.7 mM KCl.

Storage

Store at room temperature.

Applications

Prepare 1 × phosphate buffered saline.

Notes

Do not take orally.

Experiment Process

Add one bag of 1 × Phosphate Buffered Saline into 1L ddH₂O and stir to dissolve, which can be used for subsequent experiments.

*All trademarks are the property of their respective trademark owners. Some trademarks are not registered in all administrative regions.